

Actifresh Insole



Code	Actifresh - wkładki
Type	insoles
Category	Accessories
Sizes	do samodzielnego wycinania 34-46



For wide feet



Hand wash at
30 degrees



Anti-odor



Removable
insole

Description

Actifresh is made of material containing a bactericidal agent which is produced by the Swiss company, Sanitized, and which protects the foot against moulds, yeasts, algae and bacteria. The insole should be cut along the line indicated for a given size. It can be used in most shoes and boots, as well as the following models from our shop.

Additional information

SPOSÓB UŻYCIA:

1. Na tyle wkładki jest umieszczony szablon do docięcia.
Dotnij wkładkę do odpowiedniego rozmiaru.
2. Umieść wkładki w obuwiu warstwą lateksową do środka.

SPOSÓB UŁOŻENIA WKŁADKI W BUCIE:



PIELĘGNACJA:

Wkładki można prać ręcznie w maksymalnej temperaturze 30°C.

Size chart

Wkładkę należy wyciąć wzdłuż zaznaczonych linii danego rozmiaru.